

Autoclave mapping




We qualify all your static or rotary retort

°C
°F

Cold and hot spot identification


Strict compliance with acceptance criteria

Commitment* on the objective of thermal homogeneity to be achieved




Hot zone

Cold zone



*in consultation with your service quality



Autoclave mapping



We qualify all your static or rotary retort

Objective

Guarantee the industrial operation of your autoclaves in compliance with the required performance

Issue a mapping report that is also a report of approval*

Commitments

Reporting of the report within 30 days containing the applicant's requirements

Compliance of performance with acceptance criteria

Provide calibration certificates for all data loggers used

Information on the possibility of energy savings (steam, compressed air and water)

Methodology

Metrology control of the metrology of the equipment
⇒ possibility of calibration by us

Verification of controllers (temperature/ pressure)
⇒ adjustment of parameters if necessary

Critical reading of the programmed schedule
⇒ correction or if necessary re-writing our expert

Instrumentation of 5 loggers minimum per basket

*unless there is a major defect on the autoclave

